

NEW!



PIEZOSURGERY®

The original PIEZOSURGERY®
introduces now the new
implant site preparation

Implant site preparation with PIEZOSURGERY®,
the revolutionary technique – precise and efficient



- less inflammatory cells compared to implant site preparation prepared with drills
- more active neo-osteogenesis
- better control of the preparation
- less power necessary



smile. we like it!


new application for
PIEZOSURGERY

the instruments for implantology




step by step

1




OP5 CUTTING ACTION
bone perforation
CLINICAL APPLICATION
initial pilot osteotomy
(in the maxilla)
POWER Bone 1

2




IM2 CUTTING ACTION
bone perforation
CLINICAL APPLICATION
pilot osteotomy in anterior or
posterior region (in the maxilla)
POWER Bone 1

3




OT4 CUTTING ACTION
micrometric osteotomy
(about 1 mm)
CLINICAL APPLICATION
to correct pilot osteotomy axis,
to finalize the implant site pre-
paration close to the alveolar
nerve; sinus crestal approach
technique (in the maxilla)
POWER Bone 1

4



IM3 CUTTING ACTION
to enlarge the pilot osteotomy
CLINICAL APPLICATION
to finalize the implant site
preparation (in the maxilla)
POWER Bone 1

5



IMPLANT SETTING

PIEZOSURGERY®

the original

– 6 years of successful
daily use with thousands
of surgeons world-wide!

The original technology
has defined a new dimension
in bone surgery

Thanks to the controlled three-dimensional ultrasound vibrations, the original PIEZOSURGERY® technique opens up a new age for osteotomy, osteoplasty and extraction in Implantology, Periodontology, Endodontics and Surgical Orthodontics.

- Micrometric cuts maximum surgical precision and intra-operative sensitivity
- Selective cuts minimal damage to soft tissue, maximum safety for you and your patients
- Cavitation effect maximum intra-operative visibility and a blood-free surgical site



mectron Suisse
Uniline Dental Service

Le Clos, 43
1346 Les Bioux
SUISSE
tel +41 21 8454682
fax +41 21 8454714

mectron.suisse@mectron.com

mectron France
Surgytech Sarl

3706 rue de Franche Comté
39220 Bois d'Amont
FRANCE
tel +33 3 84609853
fax +33 3 84609637

mectron.france@mectron.com

mectron España
Incotrading S.A.

C° Hormigueras, 119
28031 Madrid
ESPAÑA
tel +34 91 3807490
fax +34 91 3036863
www.incotrading.net
incotrading@incotrading.net

Deutschland
Rocker & Narjes GmbH

Waltherstr. 80
51069 Köln
DEUTSCHLAND
tel +49 221 492015 0
fax +49 221 492015 29
www.rocker-narjes.com
info@rocker-narjes.com

mectron s.p.a.

via Loreto, 15/A
16042 Carasco (Ge)
ITALIA
tel +39 0185 35361
fax +39 0185 351374

www.mectron.com
mectron@mectron.com

mectron Dental (India)
PVT LTD

"Devi Plaza", #3, 9th, D' Cross
HBR Layout, Ring Road
Bangalore - 560043
Karnataka, INDIA
tel +91 80 41608671
fax +91 80 41608671
info@mectronindia.com

mectron Poland
FM Produkty Dla Stomatologii

ul. Wybickiego 1/49a
31-261 Kraków
POLAND
tel +48 12 632 19 27
fax +48 12 632 89 9
www.fmdental.pl
biuro@fmdental.pl

mectron Russia
ANAS LTD

Zemlianoi Val 2/50 KV. 31
105 064 Moscow
RUSSIA
tel +7 495 6265455
fax +7 495 917 9931
www.anas.ru
info@anas.ru

United States
Piezosurgery® Inc.

100 Matawan Road, Suite 400
Matawan, New Jersey 07747
USA
tel +1 732 335 3433
fax +1 732 335 3667
www.piezosurgery.us
info@piezosurgery.us